Español
Pulsioxímetro
8sec Rhythm

Distinguido cliente:

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English
Pulse oximeter
8sec Rhythm

Dear customer,
Many thanks for choosing to purchase a Taurus brand product. Thanks to its technology, design and operation and the fact that it exceeds the strictest quality standards, a fully satisfactory use and long product life can be assured.
Read these instructions carefully before switching on the appliance and keep them for future reference. Failure to follow and observe these instructions could lead to an accident.


Intended Use
The fingertip pulse oximeter is a portable non-invasive device intended for spot-checking of oxygen saturation of arterial hemoglobin (SpO2) and pulse rate of adult in hospitals, hospital-type facilities, and home environments.
It applies to all persons except infants and neonates.
Product is suitable for weighing more than 40kg men or women weighing more than 30kg. No special training in the use of the product, but requires the user to read the instructions before use. Readers are adults, and would be able to read and understand the operation manual; there are no reading and understanding barriers.

Operation Instruction
1. Before use, carefully read the manual.
2. Press the "power switch" key on front panel to turn the pulse oximeter power on
3. Place your one finger (forefinger or middle finger) into the cavity of the pulse oximeter. It is also available to measure if using your other fingers, such as thumb, ring finger or little finger.
4. Read the data display on the screen after displaying a stable waveform (about 3-4s). The bar graph height indicates the pulse strength.
5. The oximeter can shift the direction display by pressing the switch again. The device has two direction displays. (Fingertip pulse oximeter without direction sensor can change the direction display by pressing the switch, while the one with digital tube display cannot change the direction display.)
6. The device will turn itself off automatically after 20s idling. The finger measured has no pathological change.
Do not use the fingertip pulse oximeter in situations where alarms are required.

The device has no alarms of the SpO2 and pulse.
When the patient is lack of oxygen, we should use oximeter to analyse blood samples in order to know patient condition.
Operation of the fingertip pulse oximeter may be affected by the use of an electrosurgical unit.
This fingertip pulse oximeter is neither suitable for continuous monitoring the same site of the same patient nor suitable for measuring neonate and infant.

Check the finger touching site for every two hours. The finger measured should move or shift part if the finger skin changes. Some patients, such as the one with perfusion disorder or sensitive skin, may need frequent measurement. Long-time continuous monitoring could increase the unexpected change in skin, such as sensitive, red, blister or pressure necrosis, etc.
The results may be wrong if you did not insert the finger thoroughly in the oximeter
Please use the medical alcohol to clean the finger and the oximeter rubber touched by finger. (The oximeter rubber is medical silicon, which has no toxin, no harmful to human skin and has passed biocompatibility test.)
Do not use the device in an MRI or CT environment. It may cause errors of measurement data.
Do not use the device in the presence of flammable anaesthetic mixed with air, oxygen or nitrous oxide.
Portable and mobile RF communications equipment can affect medical electrical equipment.
These materials that contact with the patient's skin contain medical silicone and ABS plastic enclosure are all pass the ISO 10993-5 tests

for invitro cytotoxicity and ISO 10993-10 Tests for irritation and delayed-type hypersensitivity.
Temperature may rise in the process of product use, but the temperature may not exceed 40 degrees and no harm will be caused to the user.
It is suggested that the device be used indoors or in the environment where there are no strong light and too much dust.
It is best to store the product in -20°C~+55°C and ≤93% humidity.
Do not dismount, take apart or maintain the product by self. Please contact the service provider if there is any problem of the product.
If the accuracy of the measurement by the device is uncertain, check the patient's vital signs by alternative means first. Then check the oximeter. Inaccurate measurement may be caused by the following factors.
1. Ambient light radiation
2. Patient movement
3. Diagnosis test
4. Low perfusion
5. Electromagnetic, such as using the cell phone nearby
6. Electrical equipment
7. The fingernail polish
8. Artery blood is too low to measure, which is caused by shock, anemia, low-temperature or vasoconstrictor.
9. Heavy smoking patient may appear instantaneous high CO, causing the increase of the hemoglobin CO.
10. Patients with severe jaundice will have high bilirubin, which metabolizes CO that shapes significant lever carboxyhemoglobin. Thus, causes high SpO2.
When the pulse bar graph does not bounce, the displayed SpO2 or pulse rate value is potentially incorrect. Because the waveform is under-normalized, so the signal waveform when the irregular or display line, the value of SpO2 and PR displayed may be incorrect.

1. Open up the battery compartment cover
2. Install two AAA batteries into the battery compartment correctly to match the positive and negative. Slide the battery door cover horizontally along the arrow shown as the picture.
NOTE:
The device will automatically star and the luminous tube lights when the batteries installed correctly.
Please remove the batteries if the pulse oximeter will not be used for long periods of time.
Cautions:
The battery power is showed in the middle of the screen. Please replace the Batteries when battery sign shows blank space.
Caution: Batteries polarities must be installed correctly, otherwise damage might be caused to the device.
Please remove the batteries from the battery compartment if the pulse oximeter will not be used for a long period of time.
Follow local governing ordinance and recycling instructions regarding disposal of batteries, other electron devices and outer packing.
Keep the oximeter away from young children, Small items such as the battery door, battery and lanyard are choking hazards.

Maintenance
Cleaning:
Please use medical alcohol to clean the silicone touching the finger inside of oximeter with a soft cloth dampened with 70% isopropyl

alcohol. Also clean the being tested finger using alcohol before and after each test.
Do not pour or spray liquids onto the oximeter, and do not allow any liquid to enter any openings in the device. Allow the oximeter to dry thoroughly before reuse.
The device requires no routine calibration or maintenance other than replacement of batteries.
Disinfecting: Damage may be caused by disinfecting. It is recommended that the device should be disinfected when doctors think it is necessary. The device should be cleaned before disinfecting.
Maintenance: Please power up three hours per three months to prevent the damage of the water vapor, if the device does not work for a long time or is stored in a moisture environment. Additional consideration for the maintenance, please follow the hospital maintenance plan.
Clean surface of the fingertips pulse oximeter before it is used to diagnosis for patient.
Do not use the caustic or abrasive detergent. It is recommended to keep the device in dry environment. Extreme moisture may affect product lifetime and cause damage.
Do not pour and spray liquids onto the fingertips pulse oximeter, accessories, connectors, switch, or do not allow any liquid to enter openings of the oximeter to prevent the damage of the oximeter.

The use life of the device is five years when it is used for 15 measurements every day and 10 minutes per one measurement.
Please stop to use the product and follow local governing ordinance and recycling instructions regarding disposal waste.
The oximeter requires no routine calibration or maintenance during its life time. Please contact the service representative if one of the following cases occurs.

The oximeter is calibrated by simulation fingertip pulse oximeter before it leaves the factory, which can display the hemoglobin saturation.
The device cannot show normally or do not display at all.
The device cannot be powered on (New batteries have been changed.)
The sensor of the oximeter was damaged seriously and causes the oximeter stop working properly.
If the sensor is damaged (not light), indicating that the sensor has been damaged. Please do not continue to use the oximeter. And contact your local dealer for processing.

NOTE:
A functional tester cannot be used to assess the accuracy of a pulse oximeter. Clinical testing is used to establish the SpO2 accuracy. The measured arterial hemoglobin saturation value (SpO2) of the sensors is compared arterial hemoglobin oxygen (SaO2) value, determined from blood samples with a laboratory CO-oximeter. The accuracy of the sensors in comparison to the CO-oximeter samples measured over the SpO2 range of 70% - 99%. Accuracy data is calculated using the root-mean-squared (Arms value) for all subject, per ISO 80601-2-61: 2011, Medical electrical equipment— Part 2-61: Particular requirements for basic safety and essential performance of pulse oximeter equipment

A functional tester is used to measure how accurately Fingertip Pulse Oximeter is reproducing the specified calibration curve and the PR accuracy.
The model of functional tester is Index2 FLUKE simulator and the version is v3.00

Symbol definitions
⚠ - **Type BF Equipment**
⚠ - **Attention! Consult accompanying documents**
💓 - **Oxygen saturation**
💓 - **Pulse rate (BPM)**
♻ - **Symbol for disposed waste electrical and electronic equipment**

Follow local governing ordinance and recycling instructions regarding disposal of batteries
IP21 - Drip proof
CE - **CE mark**
2460
🏭 - **Manufacture's information**
🌐 - **Follow operating instructions**
EC REP - **Authorized representative in the European Community**

Factory: Shenzhen Acurio Instruments Co., Ltd.
Address : 502G Part B, Donglian Building, Chuangye 2nd Road, 23rd Subdistrict, Baoan District, Shenzhen
Phone : 86-755-21677669
EC representative: MedNet GmbH
Address: Borkstrasse 10, 48163. Münster, Germany
Tel: +49 251 32266-0 · Fax +49 251 32266-22

WARRANTY AND TECHNICAL ASSISTANCE
This product enjoys the recognition and protection of the legal guarantee in accordance with current legislation. To enforce your rights or interests you must go to any of our official technical assistance services.
You can find the closest one by accessing the following web link: <http://taurus-home.com/>
You can also request related information by contacting us.
You can download this instruction manual and its updates at <http://taurus-home.com/>

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